ORIGIN

Allodyn® is a bone allograft intended for dental surgery. The bone used is harvested from the femoral head in living donors during total hip replacement surgery. Each femoral head is harvested from a donor who has willingly consented to the use of this tissue, and is carefully selected by OST Développement prior to use.

The femoral heads thus selected are treated in our plants according to the Osteopure® process – a bone treatment process developed and patented by OST Développement

- OST Développement is a tissue bank authorised since 10 April 2001 by ANSM (the French National Agency for Medicines and Health Products Safety) (authorisation No.: FR06301T) to collect, store, transform and distribute bone allografts.

- The different variants of the Osteopure® process are authorised by ANSM under PPT Nos. 232, PPT 233, PPT 234 and PPT 235.

CONDITIONS OF SUPPLY

Allodyn® is an allograft and as such must be supplied on the basis of a prescription for a named patient use. A prescription form for a named patient use is contained inside each box. Practitioners must also inform patients of their intention to use a product of human origin.

INDICATIONS

Allodyn® is used in bone filling and augmentation grafts in post-extraction surgery, pre-implant surgery, trauma and maxillo-facial surgery.

- Socket filling
- Cyst filling
- Horizontal and vertical augmentation of the bone crest
- Sinus floor lift procedures.

STORAGE

Allodyn® does not require to be stored in a refrigerator. It can be stored at room temperature. Allodyn® does not require any special storage conditions but it must be kept in a clean, dry place, protected from direct light.

Do not use if the packaging has been opened or is damaged.

Do not use after the expiry date stated on the packaging.

Sterile, single use product. It must not be reused or resterilised.

Open the primary packaging just before use, under aseptic conditions.

TRACEABILITY

Each Allodyn® product is identified by a unique code number which links it to a single donor and a single shape. This code appears on the tracking labels inside each box containing an Allodyn® product. These labels should be attached to the recipient's records.

The implantation sheet accompanying each Allodyn® product must be completed at the time of the procedure and returned to the pre-registered address in accordance with current legislation (according to the French Public Health Code and Decree dated 14 May 2010). This feedback allows OST Développement, the original tissue bank, to fully track each Allodyn® unit.

INFECTIOUS SAFETY

OST Développement takes the greatest care to ensure the safety of each Allodyn® product. This safety is based on donor selection, the efficiency of the treatment process applied and the strict application of quality management procedures.

- Each tissue donor agrees to complete a medical selection questionnaire and to provide a blood sample which is tested for a certain

number of pathogens (HIV, HBV, HCV, HTLV and Treponema pallidum).

- The Osteopure® treatment process complies with the principles of Good Tissue Practices (Decision dated 27 October 2010) and comprises several validated steps of viral clearance / viral inactivation in particular (PASTEUR Institute TEXCELL, 1999).
- Allodyn® products are sterilised by ionising radiation (Beta rays) at a dose of 25 kGy after being placed in boxes (in accordance with ISO standard 11137).
- Each product contains tissue from only one donor and is tracked individually.
- The OST Développement quality management system is ISO 13485-certified.

CONDITIONS FOR USE

Allodyn® is rehydrated with saline solution for at least 2 minutes at the time of the procedure. Allodyn® can be used with autografts, PRF (platelet-rich fibrin) or PRP (platelet-rich plasma). The following is essential at the implantation site in order to ensure successful graft implantation:

- Good vascularisation (bleeding site)
- Close and stable contact with the graft. Osteosynthesis, which ensures the stability of corrective surgery, is required in areas of substantial mechanical constraints.
- When using Allodyn® Expand shapes, spongy or cortical-spongy blocks, pressurised hydration of the graft is required to prevent air pockets from forming within the tissue.
- Before closing the graft site, any protruding or sharp surface likely to damage the soft tissue should be removed.

CONTRE-INDICATIONS

The contraindications for Allodyn® are those that apply to grafts in general namely:

- Acute or chronic infection at the operating site
- Caution should be exercised with smokers
- Caution should be exercised with diabetic patients.

ADVERSE REACTIONS OR EVENTS

Inform OST Développement, your warehouse manager and your local Biovigilance representative immediately of any problem relating to an event or risk of an adverse reaction/event or any other anomaly observed in connection with the use of an Allodyn® product. These reports can be sent:

By e-mail to the following address:
biovigilance@ost-dev.com

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By fax: + 33 (0) 4.73.98.75.00

By telephone: + 33 (0) 4.73.98.14.28

Contact our Customer Service for more information on the Allodyn® product line

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