□ **DESCRIPTION**

OSTEOPURE is a bone allograft produced from a femoral head harvested from a living donor during hip arthroplasty. Femoral heads are validated for therapeutic use by OST DÉVELOPPEMENT, AUTHORISED Tissue Bank No. FR06301T, then processed according to the OSTEOPURE viral inactivation process so as to obtain a safe, osteoconductive cancellous allograft. The OSTEOPURE tissue preparation process applied to selected tissues has been validated for its ability to eliminate or inactivate a wide range of viruses. The products are then terminally sterilised by radiation (accelerated electron beam at a dose of 25 kGy) according to ISO standard 11137.

The tissue establishment, OST DÉVELOPPEMENT, and the OSTEOPURE tissue preparation process have been authorised by the relevant French authorities in accordance with articles 3 and 4 of Directive 2006/86/EC. The OST DÉVELOPPEMENT quality management system is ISO 13485-certified.

□ **INDICATIONS**

OSTEOPURE is used in graft procedures as a filler or as an interposition graft in orthopaedics, traumatology and maxillofacial surgery, in the following cases:

- Segmental substance loss during arthroplasty revisions
- Addition osteotomy with osteosynthesis
- Spinal arthrodesis
- Fracture with osteosynthesis
- Curettage of benign bone tumours
- Post-extraction filling
- Filling of cyst pockets
- Filling of harvested sites
- Peri-implant filling

□ **CONDITIONS FOR USE**

OSTEOPURE is supplied on personal medical prescription. The patient must be informed of the use of tissue of human origin. OSTEOPURE is a single-use product which must not be reused or resterilised. Do not use after the expiry date which is stated on the packaging.

Open the immediate packaging under aseptic conditions immediately before use; do not use if the packaging has been opened or damaged.

OSTEOPURE should be rehydrated for at least 2 minutes using physiological saline solution at the time of the procedure.

Ensure proper vascularisation of the graft site (bleeding site) and close, stable contact with the graft.

Due the cancellous nature of OSTEOPURE products, rigid osteosynthesis, per se, to ensure the stability of the surgical technique is required in areas of mechanical stress.

□ **CONTRAINDICATIONS**

Do not use in cases of acute or chronic infection at the surgical site. Do not use in weight-bearing areas without rigid osteosynthesis. Do not use for facial reconstructions.

□ **TRACEABILITY**

Each graft is identified by a unique serial number linked directly to the donor number. Traceability up to the graft recipient is ensured by the user with the traceability labels contained in each box. When required, traceability is also ensured through the implantation form accompanying each graft. This form must be completed and returned to the pre-printed address. Where applicable, traceability documents must be kept in the patient's medical record.

In accordance with current European regulations, the user must ensure that traceability records are kept for 30 years after clinical use.

□ **DONOR SELECTION**

The tissue used in the manufacture of

OSTEOPURE products is harvested from living donors who have willingly given their informed consent.

Each donor agrees to complete a questionnaire to collect information relating to his or her medical history and potentially risky practices.

In addition, each donor agrees to the testing of a sample of his/her blood, collected at the time of the procedure.



Each sample was tested for the following infectious markers:

- Anti-HIV 1 and HIV 2 antibodies
- p24 antigen
- HBs antigen
- Anti-HBc antibodies
- Anti-HCV antibodies
- Anti-HTLV I/II antibodies
- TPHA

□ **ADVERSE EVENTS OR REACTIONS**

The user must report all serious adverse events (that is, any incident related to procurement, validation, processing, storage and distribution of tissues and cells that could lead to the transmission of a communicable disease, death or risk of death, disability or incapacity in the patient, or which could result in, or prolong, hospitalisation or disease) or serious adverse reactions (that is, an unexpected reaction, including a communicable disease, in the donor or recipient related to the procurement or human use of tissues and cells that is fatal, or potentially fatal, disabling or incapacitating or resulting in or prolonging hospitalisation or a disease) within 24 hours to their OST DÉVELOPPEMENT representative.

**SYMBOLS USED:**

-  **SN** Graft identification number
-  **REF** Catalogue reference
-  **Expiry date**
-  **STERILE R** Sterilised by radiation
-  **Tissue establishment**
-  **Do not re-sterilise**
-  **Do not reuse**
-  **Storage temperature**
-  **Store in a clean and dry place, away from direct light**
-  **Do not use if the package has been opened or is damaged**

Tissue establishment:  
 OST Développement  
 13 Rue Henri MONDOR  
 63360 SAINT BEAUZIRE  
 FRANCE  
 Tel. +33 473.98.14.28  
 Fax. +33 473.98.75.00  
 service-client@ost-dev.com

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