






Dimensions

	Code	Dénomination	Long. L	Larg. I	Talon H	Pointe h	Angle	Vol. (cc)	Code LPP	Prix HT	
	OPUMCO20	OSTEOPURE COIN MDB 20/8	20	20	6	3	8,5	1,8	3317648	314,96 €	
	OPUCOCP8	OSTEOPURE COTTON 8	32	18	8	3,5	8	2,7	3338870	419,94 €	
	OPUCOPL10	OSTEOPURE LAPIDUS 10	32,4	18,6	10	4,3	10	3,6	3338870	419,94 €	
	OPUCOPL14	OSTEOPURE LAPIDUS 14	32,4	18,6	14	6	14	5	3338870	419,94 €	
	OPUCOPT10	OSTEOPURE TALUS 10			20	10	6,5	10	2,5	3338870	419,94 €
	OPUCOPT14	OSTEOPURE TALUS 14			20	14	10,5	10	3,8	3338870	419,94 €
	OPUCOCA10	OSTEOPURE CALCANEUS 10	22	23		10		4	3338870	419,94 €	
	OPUCOCA14	OSTEOPURE CALCANEUS 14	22	23		14		5,5	3324387	524,93 €	



Service Commercial :
Tél. 04.73.98.75.15 - Fax : 04.73.98.75.00
email : service-client@ost-dev.com

